

Thermal Desorption System

# TD-20




## Thermal Desorption

Thermal desorption systems trap target substances carried in mobile phase by adsorption in a sample tube packed with adsorbent. Then, by thermal desorption, the substances are delivered to a GC system. Since samples are condensed using simple operations, thermal desorption is often used when headspace methods do not provide adequate sensitivity or when gas must be collected for long periods. In addition, the sample tube is easily transported, so it is often used when samples must be collected on site.

**Major Fields of Application**

- Measurement of indoor air pollutants, air pollutants inside automobiles, clean room pollutants, and hazardous pollutants in the atmosphere
- Measurement of gases generated from parts or materials (outgassing)
- Measurement of fragrance components



## High Performance System Optimized for Microanalysis

### Simple Configuration



When a thermal desorption (TD) system is connected to a GC-MS system, the background and contamination levels must be reduced. In addition, air must be kept from leaking into the system. To satisfy these requirements, the TD-20 uses a simple flow path design and a Sulfinert secondary trap with excellent sealing characteristics.

### High Recovery Rates

The TD-20 keeps the lines from the sample tube connection port to the capillary column connection joint heated efficiently to provide uniform heat distribution. Consequently, the recovery rate of high boiling point components is high and there are almost no memory or contamination problems, even after measuring highly concentrated samples.

### Inert Lines


The TD-20 uses Sulfinert-treated material to provide inertness for all tubing in contact with gases to inhibit adsorption and decomposition. Furthermore, the TD-20 is designed specifically to couple with the GCMS-QP2010 Series, so transfer lines are short to minimize any decrease in recovery rates due to adsorption and decomposition.



## Fully Compatible with Automated Analysis

### Interlinked with Workstation Software

TD-20 functions are linked to GCsolution or GCMsolution workstation software. Not only does this allow changing conditions while serial analyses are still being performed, but it also enables analysis addition and interrupt.




GCMSolution Window for Specifying TD-20 Settings

Level#	Inj. Volu	ISTD A	Report Outp	Report File	Tuning File	Data Descri	Sampler File
1	1	1	(Level) Print				OC001.MET
2	1	1	(Level) Print				OC001.MET
3	1	1	(Level) Print				OC001.MET
4	1	1	(Level) Print				OC001.MET
5	1	1	(Level) Print				OC001.MET
6	1	1	(Level) Print				OC001.MET
7	1	1	(Level) Print				OC001.MET

For Batch Processes, Set TD-20 Parameters in the [Sampler File] Column

### 48-Sample Autosampler

The TD-20 includes an autosampler capable of holding up to 48 sample tubes. Therefore, data can be collected at night or even on weekends by batch processing in automatic operation mode.



48-Sample Sample Tray

## Exceptionally Easy to Maintain

### Parts are Easy to Replace

Maintenance tasks on the TD-20, such as replacing the O-ring seals for secondary trap or sample tubes, can be performed from the front, making maintenance a breeze. Even if transfer lines or other areas exposed to sample gases become contaminated from concentrated samples, tubing sections



can be replaced individually. In addition, joints used for the capillary column are the same as those used for the GC injection unit, making exchange or replacement easy.




Replacing the Capillary Column



Replacing the Secondary Trap Tube



Replacing O-Rings on Sealed Joints




Transfer Line Tubing

## ■ Exceptionally Easy To Operate

### Packed Cold Trap System



Before components collected in the sample tube can be introduced into the GC column, they must be focused to a band width compatible with the capillary column. Focusing can be accomplished either by the capillary cryofocus method, which cools the capillary column with liquid nitrogen or another coolant, or the packed cold trap method, which cools a trap tube packed with adsorbent. Advantages of the packed cold trap method include 1) the trap only needs to be cooled slightly lower than room temperature, which allows using electronic cooling instead of coolant, and 2) because the trap tube has a larger internal diameter than the capillary column, samples with higher moisture content can be focused without freezing the lines shut. The TD-20 uses the packed cold trap method. Therefore, it is able to obtain a more than adequate breakthrough volume, even at a relatively warm  $-20^{\circ}\text{C}$ . Consequently, electronic cooling can be used instead of time- and cost-consuming coolants, such as liquid nitrogen. Furthermore, the TD-20 is well-suited to analyzing samples with high moisture content, such as atmospheric or indoor air samples. (Tenax TA is specified as the standard adsorbent, but Carbonpak and Carbosieve adsorbents are also available for analyzing compounds with low boiling points.)



Change in Retention Times for Standard Sample Spike With 10  $\mu\text{L}$  Water

### Easy-to-Use Sample Tube


The blank level should be low when analyzing trace substances. However, the exterior of sample tubes are, for various reasons, typically contaminated with substances other than the target substance. The TD-20 sends carrier gas through the inside of the sample tube only while heating. Therefore, even if the exterior is contaminated, this minimizes the effect it can have on analysis. Also, by switching flow paths after heating is finished, the sample tube is disconnected from the carrier gas to prevent any influence from residual heat.



Effect of Applying DBP and DEHP 1000 ng to Surface of Sample Tube

### Electronic Flow Control

The TD-20 uses an electronic AFC (advanced flow controller) to control the carrier gas. Therefore, there is no need to make any adjustments to the flow rate, column inlet pressure or split ratio during thermal desorption. This also provides excellent repeatability of analytical conditions. It also allows using pressure programs or split ratio programs.



AFC-2010 Electronic Flow Controller

#### Electronic Carrier Gas Control

**AFC**

(Advanced Flow Controller)

- Does not require adjusting flow rate, column inlet pressure, or split ratio during thermal desorption
- Compatible with constant linear velocity mode
- Allows using pressure programs or split ratio programs




# TD-20 Thermal Desorption System - Analysis Examples

## Gas Samples

### Simultaneous Analysis of VOCs and SVOCs in Indoor Air


To measure indoor air pollutants, a constant-flow pump is connected downstream of the Tenax TA packed sample tube. After drawing in air for 30 minutes to 24 hours, the sample is heated for thermal desorption. Components ranging from toluene to DEHP can be analyzed simultaneously.



Upper: Indoor air from occupied studio apartment, collected for 24 hours  
Lower: 100 ng standard mixture sample

### Simultaneous Analysis of Cyclosiloxanes, Alkanes and Phthalic Esters

Cyclosiloxanes are raw materials for silicone, and trace residues are present in oil, liquid rubber and other such products. Because cyclosiloxanes are volatile, they can cause contact points in electronic parts to malfunction. Therefore, it is extremely important to control their concentration. The TD-20 is able to simultaneously measure the concentration of VOCs ranging from cyclosiloxanes to phthalic esters.






## Direct Thermal Desorption

### Comparison of Brewed Coffee and Instant Coffee


To analyze by direct thermal desorption, an empty tube is packed with a few milligrams to a few dozen milligrams of sample and placed in the TD-20 unit. The results below show the detection of pyrazines and sulfur compounds, which can affect the aroma of coffee.



## TD-20 Thermal Desorption System Specifications

Trap Tube Cooling Method	Cooled by Peltier element
Max. Number of Samples	48 tubes
Sample Tube	1/4 OD x 90 mm Tenax TA
Sample Tube Heat/Purge Flow Rate	21 mL/min to 150 mL/min (1 mL increments) at 80°C to 400°C (1°C increments)
Sampling Line Temperature	80°C to 350°C (1°C increments)
Valve Heating Temperature	80°C to 300°C (1°C increments)
Cold Trap	2 mm ID x 100 mm Tenax TA with minimum cooling temperature of 50°C below room temperature (1°C increments) If valve is heated to 260°C to 300°C, minimum cooling temperature is 45°C below room temperature and heating temperature is 80°C to 350°C (1°C increments).
Interface Temperature	80°C to 350°C (1°C increments)
Carrier Gas Control	Electronically controlled by AFC, with split ratios up to 1:200
Valve Actuation	Motor actuated
Capping/Uncapping	Controlled by stepping motor
Carrier and Purge Gases	High purity helium at 5 kPa to 900 kPa as carrier gas or dry air at 2 kPa to 300 kPa
Control Line	COM port at 9600 bps
Software	Linked TDU control software and GCMSsolution or GCsolution
Software Operating System	Windows XP, Windows Vista, Windows 7
Environment for Guaranteed Performance	Constant temperature between 18°C and 28°C (40% to 70% RH)
Power Supply	AC100 V/115 V/230 V
Size	W323 x D525 x H580 mm, not including computer
Weight	43 kg

GCMS-QP2010 Ultra + TD-20  
External Dimensions



Unit: mm

Sulfinit is a registered trademark of Silcotek Corp.  
Tenax is a registered trademark of Buchem Corp.  
Carbopak and Carbosieve are registered trademarks of Supelco Inc.



Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®".  
Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

For Research Use Only. Not for use in diagnostic procedures.  
The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.