

# Benchtop · Low Speed · Ventilated · Centrifuge

T4-L



Scan to watch product video



4000 rpm	Maximum Speed
2680 xg	Maximum RCF
Angle rotor: 6x50 ml	Maximum Capacity
≤55dB (A)	Noise Level(dBA)

## Technical Specifications:

Time range	0-99minutes adjustable	Speed accuracy	±10rpm
Power supply	AC220V,50HZ,5A,120W	Weight without rotor	18kg
Dimension	430x320x250mm(LxWxH)		

## Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	24×10/5ml	4000r/min	2680xg	PP/PC tube / Vacuum blood collection tubes
	Adapter	24×2ml			
	Adapter	24×1.5ml			
NO.2 	Angle rotor	6×50 ml	4000r/min	2680xg	PP/PC tube
	Adapter	6×15/10ml			
	Adapter	6×5ml			