

Product Description

The benchtop centrifuge MP30 combines advanced technology with high capacity.

Design

- LCD display. The centrifuge parameters can be set in operation. User-friendly program design. The adjustable modes of short spin, timing and rotational speed make the test more efficient and stable.
- Impulse centrifugation. Instant speed up to the setting value. Auto brake and fast shutdown.
- DC brushless motor with free maintenance. Faster speed acceleration and deceleration. More samples centrifugation in short time.



User friendly

- Unique appearance design with small size. Easy to use. The lid can open 90° to handle the microplates.
- Disperse the hanging liquid on tube wall quickly before or after PCR test, which improves test effects greatly.



Safety

- Produced according to international safety regulations (e.g. IEC 61010)
- Safety Electrical lid interlock, door can not open when rotor is running and machine won't run when lid open and error code appears.
- Emergency lid-lock release when power failure.